

TECHNICAL DATA SHEET

PRODUCT NAME:	Autoclave Tape without Indicator
PRODUCT CODE:	SVTWO1950, SVTWO2550
PRODUCT VARIANTS	Autoclave Tape without Indicator, 19mm x 50m Autoclave Tape without Indicator, 25mm x 50m
RELATED STANDARDS:	ISO 11140-1
INTENDED USE:	SV Autoclave Tape without Indicator is designed to stick package to be sterilized.
OPERATION CONDITIONS:	The tape should be applied to the surface of the item being sterilized before it is placed in the autoclave.
CHEMICAL COMPOSITION:	It doesn't contain any chemical component except Hot Melt Adhesive. Tape is latex free.
PRODUCT SPECIFICATIONS	SV Autoclave Tape without Indicator consists of a stickable paper roll.
STORAGE AND DISPOSAL	SV Autoclave Tape without Indicator should be stored in a cool, dry place away from direct sunlight and heat sources. It should be kept in its original packaging until ready for use. If the indicator has expired or been damaged, it should not be used, and it should be properly disposed of according to local regulations. Shelf life is 4 years.
PACKAGING	50 pcs in a box / 40 pcs in a box