

TECHNICAL DATA SHEET

PRODUCT NAME:	Crepe Paper	
PRODUCT CODE:	SVCP	
PRODUCT SIZES	Crepe Paper 300 mmx300mm	SVCP3030
	Crepe Paper 400 mm x400 mm	SVCP4040
	Crepe Paper 500 mm x500 mm	SVCP5050
	Crepe Paper 600 mm x600 mm	SVCP6060
	Crepe Paper 750 mm x750 mm	SVCP7575
	Crepe Paper 900 mm x900 mm	SVCP9090
	Crepe Paper 1000 mm x1000 mm	SVCP100100
	Crepe Paper 1200 mm x1200 mm	SVCP120120
	Crepe Paper 1500 mm x1500 mm	SVCP150150
RELATED STANDARDS:	Medical Device Regulation 2017/745	
INTENDED USE:	SV Sterilization Crepe Paper is a type of medical-grade paper that is used to wrap and protect medical equipment and supplies during the sterilization process. It is designed to withstand high temperatures and pressure during sterilization without tearing or compromising the sterility of the contents.	
OPERATION CONDITIONS:	The crepe paper is covered with autoclave tapes that allow it to be securely closed after being rolled up. Sterilization crepe paper is commonly used in hospitals, clinics, and other healthcare settings to maintain the cleanliness and sterility of medical equipment and supplies.	

PRODUCT SPECIFICATIONS

CREPE PAPER

PROPERTY	UNIT	TYPICAL VALUE
SUBSTANCE	g/m ²	60,3
THICKNESS	micron	141,00
TENSILE STRENGTH MD	N/15mm	21,3
TENSILE STRENGTH CD	N/15mm	17,1
WET STRENGTH MD	N/15mm	6,8
WET STRENGTH CD	N/15mm	5,2
BREAKING LENGTH MD	km	2,4
BREAKING LENGTH CD	km	1,8
STRETCH MD	%	9,2
BURST STRENGTH	kPa	117
AIR PERMEABILITY	ml/min	4581
WATER ABSORBENCY	g/m ²	23
pH	-	6,3
Absolute moisture content	%	5,6
STORAGE AND DISPOSAL	SV Sterilization Crepe Paper should be stored in a cool, dry place away from direct sunlight and heat sources. It should be kept in its original packaging until ready for use. If the item has expired or been damaged, it should not be used, and it should be properly disposed of according to local regulations. Shelf life is 5 years.	
PACKAGING	Packaging depends on crepe paper sizes.	